

CRYOTRON 2 - cold therapy



Cryonic cold therapy

Not-New

Information

The patented process of neuro-reflex cryo-or cold stimulation therapy works with CO₂ gas (carbon dioxide) at - 78 ° C and about 2 bar pressure, which cools the skin surface temperature of about 32 ° C within a few seconds to 2 degrees.

The CRYOTRON2

Flexible device for practice, clinic and operating room

- * Maximum mobility, stable cart with integrated gas bottle
- * Ergonomic applicator with interchangeable nozzles, laser display and infrared temperature measurement for Temperaturmessort
- * Monitor display of skin temperature, treatment time, battery status and level of the gas cylinder

Technical data:

- * Dimensions: 100 x 90 x 40 cm (HxWxD)
- * Weight: 34 kg
- * Stand-alone application by battery (12 volts)

BRAKENSIEK®
SYSTEMHAUS



Please only contact us via our registration form and only if you are a business person yourself. Every offer is subject to prior sale and includes your acceptance to our terms & conditions [Terms and conditions \(AGB\)](#). All logos and trademarks on this site are the property of their respective owners.

You will always receive an invoice including tax, which can be paid via wire transfer, in advance or cash. If you register on paypal, you may also pay by charging your paypal account with your credit card. If you are an EU member and you have a validated tax number which can be validated here: [VIES VAT number validation / MwSt.-Informationsaustauschsystem MIAS\) Validierung der MWSt.-Nummer](#) then you pay tax-free. If you are from outside the EU, you can only pay tax-free, if the export is done with the transport company, called 'Schenker'. If you are not using 'Schenker', we cannot guarantee that we will get the export confirmation in time. In that case, you have to pay tax and you get it.



**Brakensiek Systemhaus
GmbH & Co. KG
Klönnestraße 94
D-44143 Dortmund**

www.brakensiek.com

 **+49(0)231.985-00 10**

 **+49(0)231.985-2000**

BRAKENSIEK®
SYSTEMHAUS

Tim Brakensiek